



Tens Unit

Contact: Tim Grzely

Email: tgrzely@live.com

https://va-carroll.countybuyselltrade.com/health-fitness/tens-unit_madison_265244

Address:

Madison

Price:

60.00 \$

The Century 2100 is a three mode TENS unit with a timer. It can operate under three different modes. In Burst mode, the TENS unit produces a series of seven quick pulses followed by a rest. The cycle repeats itself twice every second. In Modulation mode, the TENS unit automatically varies the amplitude of the pulse every four seconds. This is a dual channel unit, which allows you to use up to four electrodes at a time (two electrodes per channel). The TENS unit is portable and may be clipped to a belt, shirt pocket, bra or other clothing.

Channels: Dual, isolated between channels

Pulse Amplitude: 0 to 80 Ma each channel

Pulse Frequency: 2 to 150Hz, adjustable

Pulse Width: 60 to 250 microseconds, adjustable

Burst Mode: 7 pulses/burst, 2 bursts/second

Timer: 30 min, 60 min or constant time

Power Source: 9V alkaline battery

Battery Life: 15 hours average usage


Size:





x



x



5



Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit

Tim Grzely
tgrzely@live.com
<https://tinyurl.com/2nok8m67>

Tens

Unit